

## Handling the slant tube/plate

Requirements for conducting these procedures in a sterile environment:

1. Store the slant tube/plate in room temperature( 25°C ) for a short-term storage after receiving.
2. Transfer the strain into a proper medium after unpacking the slant tube/plate as soon as possible ( better within a week).
3. Culture both strains from the original and the transfer tubes then incubated at a proper temperature according to the product sheet.
4. Please fill the question form and fax to BCRC if any questions including ungrown strain, contamination, culture problems and so on. **The feedback will be notified to BCRC within a month so that the strain in question can be investigated.**
5. Please feel free to contact BCRC about culture problems.
6. If you still have trouble in cultures, please assign someone to enroll the training course of "how to revive and maintain the microorganisms" hold by our institute. Please contact Mr. Tasi (886-3-5223191ext.269) for registration.